Euroline

Eurocel 01

Eurocel 01 – Your Key to a Successful Chiral Separation


KNAUER's newly developed Eurocel 01 will be a powerful tool for solving many of your chiral separation tasks. Eurocel 01 is a chiral stationary phase (CSP) based on a high quality silica matrix, coated with a derivatized cellulose polysaccharide. It is modified with the very potent chiral selector 3,5-dimethylphenylcarbamate.

The derivatization of the OH-groups results in supramolecular helical polymer structures, combining polar and π - π interactions with inclusion complexation to effect the separation of chiral compounds.

The strength of Eurocel CSPs is their impressive enantioselectivity, chromatographic efficiency and extremely high loadability. Although CSPs are typically used in normal phase mode, with Eurocel 01 it is also possible to change the chromatographic conditions to reversed phase or even 100% polar organic mobile phase mode as well. Since many chiral separations can be performed in either mode, with Eurocel 01 the user is free to select whichever eluent is most convenient.

Column: Eurocel 01,
5 μ m 250 x 4.6 mm ID
Eluent: MeOH
Flow rate: 1.0 ml/min
Temperature: 25 °C
Detection: UV

Separation of
etozolin enantiomers
using Eurocel 01 with
methanol as a polar
organic mobile phase

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Applications for Eurocel 01

CSPs have become essential for purification of chiral compounds. Polysaccharide based phases such as Eurocel 01 are the packing materials of choice for high performance chiral separations in preparative chromatography, supercritical fluid chromatography and also simulated moving bed chromatography. The main application area of CSPs is within the pharmaceutical industry.

Eurocel 01 can be applied for the purification of beta blockers as well as alkaloids, tropine compounds, amine derivatives, phenols etc.

Our applications laboratory is ready and able to assist you in designing a method for your separation task using our Eurocel 01 material.



General Properties of Eurocel CSPs

Eurocel 01 exhibits one of the highest loadabilities of all CSPs. The pressure resistance of Eurocel is comparable to that of regular silica. It is designed to be used either in normal phase mode or in reversed phase mode. Eurocel 01 is also available in a variety of particle sizes (5, 10, and 20 µm are available).

Eurocel 01 is your first choice for analytical and preparative chiral separations.

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Technical data are subject to change without notice.

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Visit www.knauer.net for details on complete HPLC systems, HPLC columns, and osmometers.