


Agilent Separation Times

2



PAGE 8 **Change Liners in Seconds!**
Innovative Flip Top System

PAGE 11 **Separate Complex LC Samples Faster**
ZORBAX Eclipse XDB-CN Column



PAGE 14 **Run Longer at Top Performance**
Preventative Maintenance Programs

CONTENTS

Events to Remember	2
ChemStation Designed for the Future	3
Find Hidden Target Compounds	4
Ask the Experts	6
Change Liners in Seconds!	8
New Life for Old LCs	10
Separate Complex Samples Faster ZORBAX Eclipse XDB-CN Column	11
Custom GC Columns for Unique Needs	13
Run Longer at Top Performance	14
Unconventional Approach for Drugs of Abuse	15





Events to Remember

**Keep Your Competitive Edge —
Stay Current with Agilent Events!**

No matter where you are in the world, Agilent makes it easier for you to keep up with the latest technology.



View a list of current events for your area!

Go to www.agilent.com/chem/events then select a country to view the local events.

**Now you can receive
Separation Times online
7 times a year - FREE!**



See a current issue and subscribe today at www.agilent.com/chem/separationtimes



KEEP YOUR LAB PRODUCTIVE

SAVE 10% WHEN YOU PURCHASE A NEW AGILENT SERVICE AGREEMENT

Your productivity depends on keeping your instruments running. When an instrument is down unexpectedly, operations come to a halt. While instruments can be repaired if a failure occurs, the time that is lost cannot be replaced. Agilent service agreements not only ensure optimum instrument uptime, they also help protect your lab's most precious asset – time. That is why Agilent Service Agreements can ultimately be the most cost-effective way to manage the maintenance and service of your Agilent instruments.

1

LIMITED-TIME OFFER: Expires January 15, 2005

2

BE PROACTIVE, NOT REACTIVE. If you have received Agilent service within the past 12 months (not under service agreement), you are eligible for this limited-time offer.

3

CALL 1-800-227-9770 IN THE U.S. AND CANADA
Or contact your local Agilent office or your participating distributor.
www.agilent.com/chem/service

**100% coverage
10% savings
0% risk**

ChemStation Designed for the Future

- Best-selling ChemStation now even better
- Single ChemStation for LC, GC and LC/MS
- Designed for future compatibility



By Chet Sharrar

GC Data Systems Product Manager

Agilent Technologies is pleased to introduce the new Revision B.01.01 ChemStation for LC, GC and LC/MS. The ChemStation Plus family is the industry's most widely sold data system, providing instrument control, data acquisition and data management. With its modular design, the system can be extended and scaled upwards as laboratory needs grow. To ensure compatibility with future PCs and their operating systems, the new revision B.01.01 uses 32-bit internal software architecture and also includes the following key improvements:

Full Support for Long File Names

Data, method and sequence file names are no longer limited to eight characters. As shown in Figure 1, long "real-world" file naming provides easier file tracking and recognition.

Additional Improvements

- Improved enhanced integrator with new integration event parameters and user interface
- Virtual GC instruments (6890 and 6850)
- Backwards compatibility for method, data and sequence files from previous ChemStation versions
- Backwards compatibility for custom macros; macro execution itself remains the same
- Supported LAN-based hardware: 1100LC, 6850 & 6890N GCs, 1100 MSD and 35900E A/D
- Supported GP-IB-based hardware: 5890, 6890A & Plus GC and 1090 LC
- Alternate 'Companion' – simplified user interface for GC production labs
- Simplified de-installation via Windows® "Add / Remove Programs" dialog
- Improved system performance and stability




Figure 1. Long file names make it easier to identify method and data files

Ordering Information

As shown in the table below, ordering Revision B.01.01 of the Agilent ChemStation for LC, GC, and LC/MS is very similar to ordering the previous revision.

Note: Existing Revision A.10.02 ChemStation products remain available concurrently with the new Revision B.01.01 ChemStation products.

Learn More

For more information, go to: www.agilent.com/chem and select DataSystems under Products & Services.

Order guide

Product Numbers:		Product Description:
Rev A.10.02	Rev B.01.01	
G2070AA	G2070BA	GC ChemStation software
G2071AA	G2071BA	Additional GC instrument control software module
G2170AA	G2170BA	HPLC 2D ChemStation software
G2171AA	G2171BA	Additional LC instrument control software module
G2710AA	G2710BA	LC/MS ChemStation software
G1656A	G1656B	ChemStation software revision for upgrading ChemStation

NOTE: It may be necessary to upgrade your PC hardware and software.

Find Hidden Target Compounds

Revolutionary MSD Software

- Save analysis time
- Increase quality of results



By Chin-Kai Meng, Ph.D. and Mike Szelewski

Typical target compound analysis requires finding target ions and meeting qualifier ion ratios. It is always very difficult to identify target compounds from high matrix background because the matrix affects the ion ratios of the target compounds. Background subtraction and manual integration are often required to be certain of the results. This is a time-consuming process; it can take an experienced analyst an hour to review/confirm one data file.

Two powerful GC/MS techniques – Retention Time Locking (RTL) and peak deconvolution – have been combined to create a quantitation and screening tool that can identify 567 pesticides and endocrine disruptors from a single run in a minute or two.

Easily Find Target Compounds in Complex Matrices

The term deconvolution is used here in the broad sense of extracting signals from a complex mixture, as depicted in Figure 1. This is the goal of AMDIS (Automatic Mass Spectral Deconvolution and Identification Software) developed by NIST (National Institute of Standards and Technology) for detecting chemicals in violation of the Chemical Weapons Convention.

Figure 2 is an example of using AMDIS to identify trace levels of one of the target compounds, methyl parathion, in grape extract. The top portion of the figure is the total ion chromatogram (TIC) of the extract. A section of the TIC (14.6 to 18.3 minutes) is expanded, as shown in the middle portion of Figure 2. It is obvious from the white TIC that it would be unlikely for an analyst to find the insecticide if the ChemStation did not report it. The bottom portion of the figure is the actual spectrum (in black) at retention time




Figure 1. Deconvolution – extraction of target signals from matrix background and interference




Figure 2. A total ion chromatogram (TIC) of grape extract and associated extracted ion chromatograms and mass spectra

(continued on page 5)

(continued from page 4)

16.57 min. The white spectral lines in the spectrum are the deconvoluted ions identified by AMDIS. Each deconvoluted/cleaned spectrum is searched against an AMDIS target library for identification. In this example, the deconvoluted spectrum matched the spectrum of methyl parathion.

Increase Confidence with Combined Report

The Deconvolution Reporting Software (DRS) for target compound analysis combines the results from Agilent ChemStation, AMDIS, and the NIST 2002 Mass Spectral Search Program (NIST02, >143,000 compounds) into one easy-to-read report, shown in Figure 3. The three result files are consolidated and sorted by both retention time and by CAS#. A printed report can be generated together with a report file in HTML format. By using RTL, detected target compounds can be disqualified when they are outside the expected retention time window. This minimizes false positives and provides a more concise final report.

The advantage of deconvolution is that it identifies target compounds in complex matrices quickly and accurately. A collaboration project with the California Department of Food and Agriculture (CDFA) Center for Analytical Chemistry demonstrated that automated DRS identified all the target compounds the analyst found in fruit, vegetable and surface water extracts. In addition, DRS found pesticides not found by the analyst. Most importantly, the time saved using the DRS was at least 30-fold, and operator concentration and skill level did not affect the consistency/accuracy of the DRS reports.

R.T.	Cas #	Compound Name	Agilent ChemStation Amount (ng)	AMDIS Match	R.T. Diff sec.	NIST Reverse Match	NIST Hit Num
4.16	95487	2-Methylphenol		87	-3.5	85	1
4.202	95487	2-Methylphenol	0.37				
4.373	108394	m-Cresol	0.23				
4.379	106445	4-Methylphenol		79	-1.1	95	1
8.754	90437	o-Phenylphenol	0.33	86	-1.0	81	1
9.966	84662	Diethyl phthalate		79	0.6	83	1
11.015	126738	Tributyl phosphate	0.33	96	5.8	92	1
16.567	298000	Methyl parathion		61	-2.6	56	1
18.414	84742	Di-n-butylphthalate	2.79	86	-0.4	94	1
19.973	10544500	Sulfur (S8)		96	-3.1	90	1
21.712	55219653	Triadimenol		74	2.3	78	1
21.948	32809168	Procymidone	2.87	97	-0.9	92	1
24.649	85509199	Flusilazole	0.34	84	4.0	85	1
26.008	563122	Ethion		95	0.4	92	1
28.413	36734197	Iprodione		88	1.7	83	1
28.635	18181801	Bromopropylate		72	1.5	85	1
29.664	117817	Bis(2-ethylhexyl)phthalate	6.49	97	1.8	91	3
30.705	13457186	Pyrazophos		78	-0.3	74	1
32.414	68359375	Cyfluthrin II		60	4.1	52	2

Figure 3. A report generated by Deconvolution Reporting Software (DRS) for the grape extract

Save Time

The deconvolution power of AMDIS increases the quality (no false positives) and confidence level (no false negatives) of the report. Compared with manual manipulation of the data by an experienced analyst (peak averaging and background subtraction of each potential target compound), the time savings using the new software program is about one hour per sample. When DRS searches for 560 target pesticide compounds, the total processing time is one to two minutes per file. Additionally, DRS has the potential to identify more compounds in complex matrices than an analyst would.

For more information

To find out more about the new Deconvolution Reporting Software, go to www.agilent.com/chem/gc. Under GC & GC/MS Data Systems, click MSD Productivity ChemStation. From there, click the link to Deconvolution Reporting Software.

To review the application note, "Comprehensive Pesticide Screening by GC/MSD Using Deconvolution Reporting Software," go to www.agilent.com/chem/library and then type publication number 5989-1157EN into the keyword field.



Ask the Experts

Information

Need technical or sales assistance? Have a technical question?

Find worldwide sales and support, online technical support or e-mail assistance. Go to: www.agilent.com/chem and select **Contact Us**. Submit your question online, directly to our experts. We'll respond promptly.



Need the best possible performance from your instruments and columns? Agilent experts have years of experience in chromatography and are happy to share their insights.

"Ask the Experts" offers practical information to make your job easier. Learn how to choose the best column for your application, optimize methods or troubleshoot problems.

How to select a guard column for an LC analysis



By **Maureen Joseph, Ph.D.**
HPLC Business Development Manager

A guard column protects the analytical column from damage caused by particulate matter and strongly adsorbed material. Ideally, the packing of the guard column should be the same as the analytical column so that the chromatography of the analytical column is not altered. If the identical material is not available, then similar or less retentive packing may be used for the guard column.

To maintain an adequate capacity for sample impurities without introducing excessive peak dispersion, choose a guard column with an internal diameter close to or the same as the analytical column and with a similar or slightly larger particle size. The guard column should be no more than 20% of the length of the analytical column. For low-volume columns, select an in-line filter for column protection.

Take the Guesswork Out of Guard Column Selection

Order your guard columns at the same time as your analytical columns. In our catalog, the LC columns and their suitable guard columns are listed together.

Or use our online compatibility guide to find which guard cartridges are compatible with each Agilent HPLC column type. For easy ordering, this guide shows pictures and part numbers of guard column hardware, with links to the online store.

Go to:

www.agilent.com/chem/columns
and select **Cartridge/Guard Cartridge Systems Compatibility Guide under LC & LC/MS Columns & Guards**.



GC/MS Tips and Tricks: Automatic Analyzer Bake-out



By Harry Prest, Ph.D.
Senior Applications Chemist GC/MS

After you clean the source and confirm that there are no leaks in the MSD, you can “bake” the source and quadrupole to rapidly lower the air-water background. With the macro programming available in the ChemStation software, you can automatically bake the MS analyzer for a period and return it to operating temperatures. You can call the bake macro in a sequence to automatically bake out the GC/MSD analyzer, return it to operating temperatures, autotune the MSD, and run a checkout sample to

confirm performance. Figure 1 shows such a sequence. Details of the macro are given in application note 5989-0678EN, available at www.agilent.com/chem/library. Type publication number 5989-0678 into the keyword field to find the application note.

Type	Vial	Sample	Method/Keyword	Data File	Comment/KeywordString
1	Keyword		Command		Macro "Bake.mac"
2	Keyword		Command		Bake 10
3	Keyword		Tune		Auto
4	Sample	1 Checkout sample	CHECKOUT	Checkout1	test of system performance
5	Sample	1 Checkout sample	CHECKOUT	Checkout2	test of system performance-rep
6					
7					
8					
9					
10					
11					
12					

Figure 1. Sequence to bake analyzer, tune and confirm performance

Choosing between on-column and splitless GC injection



By Eberhardt Kuhn, Ph.D.
Senior Applications Chemist

Both on-column and splitless injection techniques are suitable for analyzing trace analytes (analytes in the part-per-million level or less in the final injected sample). Each technique offers advantages over the other, and frequently the choice of technique is a compromise between a variety of variables.

Regardless of injection technique, regular maintenance of the column and injection port liner is very important. Even “clean” samples, such as water extracts and pure solvents, contain some nonvolatile material.

This contamination, as it accumulates from numerous sample injections, may cause chromatography problems such as activity or retention time shifts. Regular maintenance prevents these problems.

Which injection technique is best for your method?

Variable	Common Examples	Yes	No	Preferred Technique
Sample residues	soil/biological extracts, buffers, sample degradation products, salts	x		Splitless
			x	Both
Active analytes	-NH _x , -OH, -COOH, -CHO, -SH, -NHOH, cyclic -N=.	x		On-column
			x	Both
Thermally labile analytes	epoxides, ureas, explosives	x		On-column
			x	Both
Low volatility analytes	motor oil, waxes, triglycerides	x		On-column
			x	Both



Information

To select the right liner, visit our Web site at <http://www.chem.agilent.com/cag/cabu/linerselect.htm>

For an instructional video on how to properly install a liner, go to <http://www.chem.agilent.com/cag/cabu/howtovideo.htm>



Change Liners in 30 Seconds!

Innovative Flip Top System



Eliminate the time and dangers of handling heated parts; avoid frustrating searches for special wrenches or tools; decrease downtime; and extend column life by minimizing exposure to ambient air.

Exclusively from Agilent, the Flip Top Inlet Sealing System is great news for anyone who changes inlet liners on an Agilent 6890, 6850 or 5890 GC or GC/MS.

The Flip Top Inlet Sealing System allows you to safely and reliably:

- **Change an inlet liner in as little as 30 seconds**
- **Without any tools**
- **Without any leaks**
- **No handling of heated parts, which means no burns**
- **No lengthy downtime, for increased productivity**
- **Minimal exposure to ambient air, for longer column life**

Simplifying a Frequent Maintenance Task

When running a GC or GC/MS system, a hot inlet may have to be opened as often as every 24 hours to change a dirty inlet liner. The special wrench used for this operation is often misplaced, bent at an odd angle, too thin, or simply awkward to use. Even after cooling, the inlet nut is usually too hot to handle so it must be turned several times with the wrench before it will release the top weldment assembly of the injection port.

Once the dirty liner is changed, the inlet nut must be replaced on the injection port and occasionally the wrench slips off the nut, leading to scrapes, burns and cuts. It can take five minutes or more just to replace the liner. Finishing the procedure, including re-equilibration, can result in **15 minutes of total downtime.**

Changing liners on a GC/MS causes even more problems as ambient air is drawn into the capillary column, through the hot MS interface, and into the heated source. The result: multiple problems for the operator, including shortened column life and air background in the MS.

Order guide

Description	Part Number
Flip Top Inlet Sealing System	5188-2717
Inlet Liner O-Rings for the Flip Top System, 10/pack	5188-2741

How It Works

- A levered arm attaches to any Agilent 6890, 6850 or 5890 GC insert weldment and locks to the injection port using an adapter ring screwed onto the inlet.
- You simply lift the arm of the Flip Top, which releases the insert weldment from the injection port and allows instant access to the liner.
- You reverse the process to re-seal the weldment to the port.

You've Got to See It to Believe It

How easy is the new Flip Top Inlet Sealing System to use? Just view our online video demonstration of the Flip Top in action at:

www.agilent.com/chem/fliptop4u

With its inherent safety, simplicity and reliability, the Flip Top Inlet Sealing System is perfect for anyone who changes liners on an Agilent GC.

OTHER INLET ADVANCES:

New GC Inlet Septa Packaging for Greater Cleanliness and Savings

Innovative new tri-fold blister pack guarantees that the last septum is as clean as the first:

- Individually sealed septa for ultimate cleanliness.
- No worry about contaminating other septa in a common container.
- Easily dispensed one at a time.
- No septa clumping or sticking.
- Easy to see when to reorder.
- Package fits easily into drawers.
- Convenient quantities of 50 or 100.
- High-quality PET packaging — tested using GC/FID, GC/MS and GC/ECD to ensure no interfering background peaks.

New economical packs save you money:

- Economical 50/ or 100/pack septa replace all of the general purpose 24/ and 25/pack septa.
- 50/pack septa offers savings from 10% to 25% over the old 24/ or 25/pack septa. You save even more with 100/pack septa!

Order online

For price and quantity listing, visit www.agilent.com/chem/store

or call your local Agilent sales office or authorized Agilent distributor.

1-800-227-9770



New Life for Old LCs

Factory Refresh Service

New service! We're constantly improving the quality of Agilent parts. When you have your 1100 Series LC system refreshed, you get the latest, most productive technology.

Even your rugged, reliable Agilent 1100 Series liquid chromatograph — the world's top-selling LC system — can benefit from an occasional factory refresh. Our 1100 Factory Refresh Service will return a system that has performed for years to current factory specifications and performance, allowing complete confidence in your analyses and enhanced efficiency.

Agilent's 1100 Factory Refresh Service includes:

- Replacement of all consumable parts in the solvent path
- Preventive maintenance
- Lamp replacement
- Firmware update
- Everything necessary to bring your 1100 LC up to current new-instrument shipping specifications

Who better to refresh your 1100 LC system than the original manufacturer? You can enjoy peace of mind knowing that our factory-trained technicians use state-of-the-art, traceable diagnostic equipment to ensure that your 1100 LC is running to factory specifications.

This service is available only in Europe and North America. Find out more by calling 1-800-227-9770 in the U.S. or Canada. In Europe, please contact your local Agilent sales office. Or go online to learn more about Agilent's service offerings and special promotions at www.agilent.com/chem/service.

Why settle for ordinary service when you can receive the service intended for your extraordinary 1100 LC system?

Which Service is Right for You?

Service Description	Agilent 1100 Factory Refresh Service	Non-Agilent Service Provider
Up-to-date factory-trained service technicians	√	NO
Service performed to factory specifications using traceable diagnostic equipment	√	NO
Testing to factory specifications	√	NO
Replacement of detector lamp	√	MAYBE
Firmware upgrade to latest revision	√	NO



Separate Complex LC Samples Faster

ZORBAX Eclipse XDB-CN Column

- Use ZORBAX Eclipse XDB-CN for polar and non-polar analytes
- High resolution, efficient separations



By Maureen Joseph, Ph.D.
HPLC Business Development Manager

It can be difficult to attain the ideal separation of a complex sample, such as a mixture of pesticides. This type of mixture often contains both polar and non-polar compounds. An HPLC separation on a C18 reversed phase column, such as the Eclipse XDB-C18 column, is frequently long and inefficient. For this type of sample, a more polar bonded phase is a powerful choice. A more polar, cyano bonded phase, such as the new Eclipse XDB-CN, can reduce the retention times of the most non-polar analytes while still maintaining good retention of the more polar analytes.

Figure 1 shows the separation of a mixture of six pesticides on the Eclipse XDB-C18 and the new Eclipse XDB-CN column. On the Eclipse XDB-C18 column (Fig. 1A), all peaks are well resolved, but the analysis time is 54 minutes due to the excessive retention of the pesticide pencycuron. When the pesticides are separated using the more polar selectivity of the new Eclipse XDB-CN column (Fig. 1B), **the analysis time is reduced by 85% under the same mobile phase conditions**. This is primarily because the retention of the most non-polar pesticide is reduced. At the same time, high resolution of all components and the trace impurities has been maintained.


Figure 1C shows that while the retention time of the last peak can be reduced on the Eclipse XDB-C18 to match that of the Eclipse XDB-CN column by increasing the organic in the mobile phase, the *k* values indicate that all other peaks are less well retained and the resolution is inferior. Therefore the new Eclipse XDB-CN column, with its more polar bonded phase, provides the best separation of this mixture of polar and non-polar pesticides.

For more information


Download free application notes about the Eclipse XDB-CN column! Go to www.agilent.com/chem/library.

At the top of the page, enter the keyword Eclipse XDB-CN and click the Search button.

Urea Pesticide Analysis Comparison of Eclipse XDB-CN vs. Eclipse XDB-C18



Peak	Capacity Factors (<i>k</i>)		
	A	B	C
1. Fenuron	0.71	0.53	0.31
2. Monuron	1.75	0.90	0.56
3. Monolinuron	2.65	1.19	0.74
4. Diuron	1.00	4.19	1.49
5. Linuron	6.29	1.97	1.34
6. Pencycuron	36.1	3.76	3.79





Two Great Discounts and Two Great Financing Options Make it Easy to Get the Latest Technology

Agilent instruments have always been a sound investment. Now we give you the Agilent Value Promise of at least 10 years of guaranteed value from your Agilent instruments ... plus the added value of significant discounts.

TWO GREAT DISCOUNTS:

1

Trade up to new Agilent technology and save 15%. If you own an Agilent 5890 GC, 5970 GC/MSD, 5971 GC/MSD, 5972 GC/MSD, 1050 LC, 1090 LC, or 4500 ICP-MS, you get a 15% discount when you buy or lease one or more Agilent 6890 GC, 6850 GC, 5973 inert MSD, 1100 Series LC, or 7500 ICP-MS systems. For details, visit www.agilent.com/chem/library, then enter 5989-0640EN in the keyword field.

2

Buy or lease any 6 new Agilent systems and get the average price of one system back. If you own an Agilent 5890 GC, 5970 GC/MSD, 5971 GC/MSD, 5972 GC/MSD, 1050 LC, 1090 LC, or 4500 ICP-MS, you get a substantial 17% discount when you buy or lease 6 or more 6890N GC, 6850N GC, 5973 inert GC/MSD, 1100 Series LC, or 7500 ICP-MS systems. All support options, sampling systems and data systems are included. For details, visit www.agilent.com/chem/library, then enter 5989-0639EN in the keyword field.

TWO GREAT FINANCING OPTIONS:

1

OPTION 1: Lease to Own. You pay a small amount monthly or quarterly, then take title when you make the last payment. You can bundle hardware, maintenance and services into one total, economical package, with no upfront capital expense.

2

OPTION 2: Technology Refresh. You can lease Agilent products for various periods of time, paying less per month than with a traditional lease. At the end of the contract, you can return the equipment, continue to lease month-to-month, upgrade to new equipment with a new lease, or buy your leased equipment at fair market value. This lets you stay current with the latest technology without the long-term costs of equipment ownership. For details, visit www.agilent.com/chem/library, then enter 5989-0641EN in the keyword field.

8 Excellent Reasons to Move Up Within the Agilent Product Line

- The Agilent name is your assurance of unequalled dependability, reliability, ease of use and enhanced productivity.
- Optimize your laboratory's performance.
- Maximize your return on investment.
- Unrivaled technical support and educational resources to improve your laboratory's efficiency.
- A comprehensive array of instruments, supplies and services, each rigorously designed, tested and manufactured to the industry's highest standards.
- Improve analytical methods while maintaining compatibility and support.
- Enjoy familiar instrument operation and service, with minimal risks involved with change.
- Keep common instrument supplies for easier support and budget management.

A PerfectFit Bonus

With any of these offers, you'll also receive Agilent's PerfectFit 15% discount on an unlimited number of consumables and supplies when you make your purchase.

Custom GC Columns for Unique Needs

For the rare times when you need a capillary GC column that is not among the nearly 1,800 regularly available from Agilent, our Custom Column Shop can usually make the column for you. Some columns can be on their way in only a few days while others (such as custom columns with unique film thicknesses or diameters) may take up to three to four weeks for manufacturing and testing. Deactivated fused silica tubing and CE columns are also available from the Custom Column Shop.

New Web site for Custom Column Quote

NOW to speed the process, customers and distributors in the U.S., Canada and Puerto Rico can submit requests for a custom GC column quotation through our Web site. Go to:

www.agilent.com/chem/gc.customcolumn

After you submit your requirements and request a quotation, within two business days an Agilent Call Center representative will either fax or e-mail a quotation to you. The representative normally will follow up with a telephone call to make sure that the quotation arrived, to answer any questions, and to accept an order if you want to proceed.

Customers and distributors outside of the U.S., Canada and Puerto Rico should continue to contact their local Agilent sales office to submit a request for a custom column quotation.

A convenient order form is available on page 747 of the 2002-2003 Chromatography and Spectroscopy Supplies Catalog (5988-4785EN).

Information


Unsure what custom column you need or whether you need one at all?

Contact Technical Support or your local Agilent Technologies representative.
Go to: www.agilent.com/chem and select **Contact Us**.



Custom GC Column Types

	Column type
√	Individual custom column
√	Individual column with fused silica guard column or transfer line (including DuraGuard integrated guard columns)
√	Dual column assembly (two columns attached end to end or in parallel with a Y connector)
√	Standard column rewound on a 5" or custom cage
√	Custom fused silica tubing (no stationary phase)
√	Other



Run Longer at Top Performance

Preventive Maintenance Programs

- Ensure optimum performance
- Avoid unplanned downtime
- Extend instrument life



By Gary Lee


Services Marketing Specialist

In today's laboratory environment, technical professionals require their instruments to provide highly sensitive and accurate analyses reliably — day in and day out. Important factors that can have an impact on instrument performance, quality of results and instrument lifetime include routine care and regularly scheduled preventive maintenance (PM). Performing regular preventive maintenance will help to assure optimum performance and avoid expensive unplanned downtime and repairs, as well as substantially extend the useful life of your equipment. This is especially true for:

- Method development and R&D labs where instruments are subject to a diversity of analyses and, sometimes, instrument users.
- Analytical services labs where instruments need to be available to analyze a broad range of samples.
- Quality assurance labs that require accurate analytical data on a timely basis. In many cases, analytical systems in this environment run on a 24-hour basis and must operate within stringent regulatory guidelines. Instrument to instrument uniformity and performance are critical.

Regularly Scheduled Preventive Maintenance reduces downtime and repair cost.

Studies of instrument repair call activities conclude that over 60% of repairs are not functional failures but are mainly due to lack of preventive maintenance — replacing components that wear on a regular basis. A regularly scheduled preventive maintenance program is a predictable expenditure that reduces unplanned instrument downtime and repair cost, as well as inconvenience and lost productivity. It increases the effectiveness of laboratory personnel.



The graph above shows what happens as a lab grows with new equipment. Failure rate increases as the instrument population increases because most modern instrument failures occur in the first year or two of instrument life. As the population stabilizes, the failure rate should flatten (predicted failure rate extrapolation). However, when the laboratory implements a thorough maintenance program, the failure rate decreases by up to 25% for mechanical systems.

Developing a preventive maintenance program

The frequency of preventive maintenance largely depends on your application requirements, workload and regulatory requirements. In developing a preventive maintenance program for your instruments, consider the following factors:

- Type of instrument — Gas and liquid chromatographs, instruments with liquid pumping systems, auto-samplers and mass selective detectors etc. generally require more preventive maintenance than simpler instrument systems such as UV/VIS or IR.
- Run-time requirements — This can be based on hours used or number of injections (e.g. 24-hour operations). Preventive maintenance should be performed at least twice a year for QA/QC labs to maximize uptime. In FDA/GLP environments, it is general practice to perform preventive maintenance prior to operational qualification.

- Instrument usage — If your analytical instrument is used to analyze a wide range of samples, uses several analytical conditions and is run by multiple users, then regular preventive maintenance minimizes downtime and time spent troubleshooting instrument problems. This helps free scientists to develop methods, and to focus on R&D and/or analyzing samples. Two PM calls per year may be required depending on the samples and operating conditions.
- Instrument age — As analytical instruments age, regular preventive maintenance will keep them operating to manufactured specifications.

Qualified Personnel Should Perform Preventive Maintenance using Genuine Parts

Instrument damage can result when unqualified personnel attempt to provide preventive

(continued on back page)

Unconventional Approach for Drugs of Abuse

- Novel approach uses LC/MS instead of immunoassay and GC/MS
- Faster sample turnaround
- Flexibility for new analytes



By Ron Majors, Ph.D.
Senior Applications Chemist

Located just below the Arctic Circle in Norway, the city of Trondheim seems an unlikely place to find one of the largest drugs of abuse (DOA) testing labs in Europe, the St. Olavs Hospital/Trondheim University Hospital. Under the auspices of Kolbjorn Zahlsen of the Department of Clinical Pharmacology, the laboratory has adopted an analysis strategy unlike that of other drugs of abuse testing laboratories throughout the world. Most laboratories do a preliminary screen of urine samples using immunological methods. If they find a positive result for drugs, they retest the urine using gas chromatography/mass spectrometry (GC/MS) for confirmation. Zahlsen's approach is to perform screening using only liquid chromatography/mass spectrometry (LC/MS), with specificity and performance that is comparable to GC/MS confirmation.

Cost-Effective Analyses

It might seem that LC/MS would not be cost-effective in handling 3,000 analyses per day in samples from all over Norway. However, the LC/MS technique competes effectively with the expensive immunochemicals and the rigid testing protocols of the immunoassay-GC/MS approach.

Immediate Answers without Rework

There are other reasons to prefer LC/MS as the analysis technique. Although immunoassay is fast and simple, there are inherent difficulties with the technique. Most of the tests are not specific, broadly looking at DOA classes such as amphetamines, opiates, benzodiazepines, barbiturates and so on. If a positive result is obtained, the urine sample must be reanalyzed; the suspected analytes are chemically derivatized with a reagent required for GC volatility, and then the reanalysis is performed by GC/MS. With LC/MS, the specific compound is determined in the initial analysis without the need for further rework. Such specificity is especially valuable when a urine sample con-

tains several classes of illegal drugs in combination with legal drugs, and has proven very useful for DOA control of samples from methadone treatment maintenance. For some categories of samples, such as forensic samples and samples that have serious legal implications, GC/MS confirmation is performed in addition to the LC/MS screen.

Short Turnaround Time

Another advantage of Zahlsen's strategy is short turnaround time, crucial in drug overdose situations. In the clinical toxicology department, typical turnaround times are one to two hours after the sample is received — which sometimes can mean life or death for the affected individual. Sometimes, GC derivatizations alone in the conventional DOA laboratory can take an hour. Then additional time is required for sample preparation, immunoassay and the GC/MS analysis.

Fast Method Development

LC/MS is superior for analysis of new drugs. Development of the antibodies for a new immunoassay may take up to a year. With LC/MS, the method can be developed in one or two days if a reference substance is available.

It is this flexibility that has allowed Zahlsen's laboratory to extend its LC/MS work to ethical drugs (pharmaceuticals) in the therapeutic drug-monitoring (TDM) program. The ability to monitor a patient's drug levels in a hospital setting is crucial to maintain quality medical treatment. Patients often require different optimum levels of drugs so that they do not suffer from adverse effects of having too much or too little active drug in their systems. Over 70 drugs are now routinely monitored in patient samples.

Success Leads to Laboratory Growth

Zahlsen started method development with mass spectrometry in 1990 with one

Hewlett-Packard GC/MS system. He bought his first LC/MS system, an 1100 Series LC/MSD system, in 1998. At that time he had 15 employees handling about 50,000 analyses in urine and serum samples. Last year, as depicted in Figure 1, the laboratory performed




Figure 1.

820,000 analyses with 80 employees and 10 Agilent GC/MS systems, 24 Agilent LC/MSD systems and one Agilent LC/MSD Trap. This year, he expects to perform more than a million LC/MS analyses in his lab. Kolbjorn is shown below with one of his LC/MSD systems.



Kolbjorn Zahlsen can be contacted at the email address: zahlsen@online.no.

Run Longer at Top Performance

(continued from page 14)

maintenance on analytical instruments. Improper installation or non-compatible consumables used for preventive maintenance can extend downtime of your instrument. Agilent engineers have extensive training and many years of experience in identifying potential problems. They can propose appropriate actions that need to be taken to optimize performance and avoid system downtime.

Agilent has specifically engineered its line of instruments with detailed, consistent preventive maintenance protocols to assure fast and

accurate alignment and adjustments. You need more than just a tune-up or a general inspection of your instruments. Preventive maintenance specifically designed for your instrument ensures that you will be able to take appropriate action before any potential downtime occurs. Agilent's PM protocols specify genuine Agilent replacement parts and consumables to assure optimum performance.

Select an Agilent PM program to meet your needs.

Agilent offers a variety of preventive maintenance programs performed by highly trained

and certified service engineers to meet your unique requirements. You can configure your Agilent preventive maintenance program with a comprehensive service agreement, industry specific service bundle, or a PM-only agreement. Or purchase PM on a per-incident basis. At Agilent, we use factory-engineered checklists of preventive maintenance procedures for consistent service from laboratory to laboratory. Minimize your unplanned downtime by trusting Agilent to handle your preventive maintenance needs. Call today! For more information, see www.agilent.com/chem/service.

Save Money on Instruments, Supplies, and Services — Take Advantage of Special Offers from Agilent!

Intensifying workload? More demanding analyses?

To meet today's challenges, you need the best analytical tools available. Our limited-time special offers give you the opportunity to enjoy the quality and dependability of Agilent products at the lowest possible prices.

Benefit from the latest special offers from Agilent!

Visit www.agilent.com/chem. Under **Additional Information**, click **Special Offers**.



Agilent Technologies

AGILENT TECHNOLOGIES
LIFE SCIENCES AND CHEMICAL ANALYSIS GROUP
71 SOUTHGATE BLVD
NEW CASTLE, DE 19720-9939

PRESORTED
STANDARD
U.S. POSTAGE
PAID
WILMINGTON DE
PERMIT NO 1387

This information is subject to change without notice.

© Agilent Technologies, Inc. 2004

June 15, 2004
5989-1198EN