

For multi-mycotoxin analysis

18 mycotoxins at one go!

The CrossTOX® columns by LCTech allow a high efficient sample clean-up of regulated and expected mycotoxins in many matrices. At the same time, the columns improve the conventional dilute-and-shoot process via an QuEchERS based procedure:



High purity of extracts due to special sorbents for high matrix binding; no separate syringe filtration



Significantly less 13C Internal standard required due to special affinity mechanism of the columns



Extension of cleaning and maintenance intervals of your LC-MS/MS system due to the high purity of the extracts



Only one column for a wide range of applications

Catch them all

Aflatoxine B1, B2, G1, G2

Ochratoxin A

Zearalenon

Deoxynivalenol

Fumonisine B1, B2

T-2

HT-2

Nivalenol

3-Acetyl-DON

15-Acetyl-DON

DON-3GIc

Sterigmatocystin

Citrinin

Diacetoxyscirpenol



The CrossTOX® columns are suitable for both manual and automated processing. Please find in the following your advantages compared to the conventional Dilute-and-Shoot procedure:



LCTech






Range of Applications

The range of applications for CrossTOX® is remarkably wide. With only one extraction protocol the columns can be used for all matrices commonly in this analysis, e.g. cereals, nuts or dried fruits. The UV-chromatogram shows the very high clean-up efficiency of CrossTOX®. The **pink** line represents a wheat sample cleaned via the conventional dilute-and-shoot procedure (syringe filtration). The **black** line represents clean-up via CrossTOX®.

This shows the substantial improvement using dilute-and-shoot compared to CrossTOX® column clean-up.

*A constantly growing list of all matrices tested so far can be found at www.LCTech.de. We would also be happy to receive your feedback.

Recovery rates compliant to EN 401/2006




Processing via FREESTYLE

For automated processing of the CrossTOX® columns the FREESTYLE system offers two options:

- from raw extract to cleaned sample or
- from raw extract to chromatogramm (full automation)

In order to realize the highest possible degree of automation, the FREESTYLE SPE robotic systm is simply upgraded with an HPLC Direct-Injection module for direct injection into your LC-MS/MS system. This enables the parallelization of clean-up and analysis for up to 120 samples/day and thus to an optimum utilization of your LC-MS/MS system.





Ordering Information

P/N	Description	Pc. / Pck.
17900	CrossTOX® Clean-up Columns Manual processing	100 Pc. / Pck.
17901	CrossTOX® Clean-up Columns Automated processing	100 Pc. / Pck.
FREE -STYLE Robotic System FREESTYLE SPE + HPLC Direct-Injection-Module. We are happy to send you a suitable configuration on request.		

