


NEW

YMC
LF-0023E

Triart C18 ExRS

UHPLC & HPLC Columns

Particle size
1.9 μm
3 μm
5 μm



YMC CO., LTD.
<http://www.ymc.co.jp>



**CHROMATOGRAPHIC
SPECIALTIES INC.**
www.chromspec.com

1-800-267-8103 • sales@chromspec.com • tech@chromspec.com



Triart C18 ExRS

UHPLC & HPLC Columns

Alternative selectivity to standard C18 columns

Features

- C18 phase bonded with high density on inorganic/organic hybrid silica gel
- Excellent selectivity of isomers and structural analogs
- Superior chemical durability
- Identical selectivity over particle sizes


Specifications

Base	: Organic/inorganic hybrid silica
Particle size	: 1.9, 3, 5 μm
Pore size	: 8 nm
Specific surface area	: 430 m^2/g
Bonding	: Trifunctional
Carbon content*	: 25%
Endcapping	: Yes
pH range	: 1-12


*Containing 8% for hybrid silica base material

Alternative selectivity to standard C18 columns

[Comparison of hydrophobicity and hydrogen-bonding capacity of various columns]





Comparison of fundamental separation selectivity




A mixture that consists of compounds with various characteristics is analyzed with Triart C18 and Triart C18 ExRS. Triart C18 ExRS has lower polarity and higher hydrophobicity than the standard C18 column, Triart C18. It also shows improved planar cognitive ability.

Ideal for separations of structural analogs




Column : 5 μ m, 150 X 3.0 mml.D.
 Eluent : THF/acetonitrile (10/90)
 Flow rate : 0.425 mL/min
 Temperature : 30°C
 Detection : UV at 265 nm
 Injection : 4.25 μ L (10 μ g/mL)



Triart ExRS is effective for separation of structural analogs. This feature is especially preferable for separation of pharmaceuticals and structurally similar impurities.

Improved durability



Column : 5 μ m, 150 X 4.6 mml.D.
 Eluent : 50 mM K_2HPO_4 - K_3PO_4 (pH 11.5)/methanol (90/10)
 Flow rate : 1.0 mL/min
 Temperature : 40°C
 Sample : benzyl alcohol

High density bonding of C18 greatly contributes to improvement of chemical durability.


Ordering information

Particle size (μm)	Pore size (nm)	Column size inner diameter X length (mm)	Product number
1.9	8	2.0 X 20	TAR08SP9-0202PT
		2.0 X 30	TAR08SP9-0302PT
		2.0 X 50	TAR08SP9-0502PT
		2.0 X 75	TAR08SP9-L502PT
		2.0 X 100	TAR08SP9-1002PT
		2.0 X 150	TAR08SP9-1502PT
		3.0 X 50	TAR08SP9-0503PT
		3.0 X 75	TAR08SP9-L503PT
		3.0 X 100	TAR08SP9-1003PT
		3.0 X 150	TAR08SP9-1503PT

Particle size (μm)	Pore size (nm)	Column size inner diameter X length (mm)	Product number
3	8	2.1 X 20	TAR08S03-02Q1PTH
		2.1 X 33	TAR08S03-H3Q1PTH
		2.1 X 50	TAR08S03-05Q1PTH
		2.1 X 75	TAR08S03-L5Q1PTH
		2.1 X 100	TAR08S03-10Q1PTH
		2.1 X 150	TAR08S03-15Q1PTH
		3.0 X 50	TAR08S03-0503PTH
		3.0 X 75	TAR08S03-L503PTH
		3.0 X 100	TAR08S03-1003PTH
		3.0 X 150	TAR08S03-1503PTH
		4.6 X 33	TAR08S03-H346PTH
		4.6 X 50	TAR08S03-0546PTH
		4.6 X 75	TAR08S03-L546PTH
		4.6 X 100	TAR08S03-1046PTH
		4.6 X 150	TAR08S03-1546PTH
		4.6 X 250	TAR08S03-2546PTH
5	8	2.1 X 20	TAR08S05-02Q1PTH
		2.1 X 33	TAR08S05-H3Q1PTH
		2.1 X 50	TAR08S05-05Q1PTH
		2.1 X 75	TAR08S05-L5Q1PTH
		2.1 X 100	TAR08S05-10Q1PTH
		2.1 X 150	TAR08S05-15Q1PTH
		3.0 X 50	TAR08S05-0503PTH
		3.0 X 75	TAR08S05-L503PTH
		3.0 X 100	TAR08S05-1003PTH
		3.0 X 150	TAR08S05-1503PTH
		4.0 X 150	TAR08S05-1504PTH
		4.0 X 250	TAR08S05-2504PTH
		4.6 X 33	TAR08S05-H346PTH
		4.6 X 50	TAR08S05-0546PTH
		4.6 X 75	TAR08S05-L546PTH
		4.6 X 100	TAR08S05-1046PTH
4.6 X 150	TAR08S05-1546PTH		
4.6 X 250	TAR08S05-2546PTH		

※Please inquire us for a product other than listed above.

Consideration of connector and column fittings



The above information is as of October, 2014 and is subject to change without notice. ©2013 YMC CO.,LTD.All Rights Reserved.

141027A



**CHROMATOGRAPHIC
SPECIALTIES INC.**
www.chromspec.com

1-800-267-8103 • sales@chromspec.com • tech@chromspec.com

