

The **EXP2** Pre-Column and In-Line Hand-Tight Filters are ideal for protecting small-particle UHPLC columns and systems at extreme pressures.

**EXP2** Filters help extend the life of a UHPLC system and protect your investment without sacrificing performance.

**EXP2** Filter Kits are more than 50% off MSRP



Pre-Column Filter Holder



In-line Filter Holder



EXP2 Filter Cartridges

## Advantages of **EXP2** Filters

- Provides a better, more robust, and cost-effective product compared to alternatives
- Extends the life of UHPLC system and protects your investment without sacrificing performance
- Holders are reusable and ensure auto-adjusting zero-dead-volume connections for all 10-32 ports
- Low-volume, low-dispersion cartridges are uni-directional
- Free-Turn® architecture allows the user to change cartridges by hand without breaking fluid connections on the holder inlet/ outlet.
- Hardened stainless steel end cap eliminates galling
- All holders include EXP Fittings
- Approved for use at 20,000+ psi (1,400+ bar)


## EXP2 FILTER KIT PROMO DETAILS

PART NUMBER	KIT CONFIGURATIONS
15-04-05497	EXP2 Pre-Column Holder and (3) Filter Cartridges, 2µm
15-04-05498	EXP2 In-Line Filter Holder and (3) Filter Cartridges, 2µm
15-04-05499	EXP2 Pre-Column Holder and (3) Filter Cartridges, 0.5µm
15-04-05500	EXP2 In-Line Filter Holder and (3) Filter Cartridges, 0.5µm
15-04-05501	EXP2 Pre-Column Holder and (3) Filter Cartridges, 0.2µm
15-04-05502	EXP2 In-Line Filter Holder and (3) Filter Cartridges, 0.2µm



**OPTIMIZE** TECHNOLOGIES

www.optimizetech.com



**CHROMATOGRAPHIC  
SPECIALTIES INC.**

www.chromspec.com

