

CHROMSPEC™ SYRINGE FILTERS

CHROMSPEC™ . . . The Name for Affordable Quality

- Ultra clean
- Autoclavable
- Fusion sealing ensures purity of filtrate and eliminates bursting
- Best combination of large surface area with minimal hold-up volume
- Durable polypropylene housing for increased solvent resistance
- Female Luer Lock inlet connection to ensure a secure fitting and a Male Slip Luer Fitting
- Individually labeled to prevent confusion or unnecessary waste
- Available in 4 mm, 13 mm and 25 mm sizes for an increased selection for your filtering needs
- Samples available



CHROMSPEC™ NON-STERILE SYRINGE FILTERS

| Size (mm) | Membrane Type | Pore (µm) | Product Number | Pkg. 100 Price | Product Number | Pkg. 500 Price | Product Number | Pkg. 1000 Price |
|-----------|---------------|-----------|----------------|----------------|----------------|----------------|----------------|-----------------|
| 25 | Nylon | 0.20 | CS25N201 | \$140.00 | CS25N202 | \$635.00 | CS25N203 | \$1,050.00 |
| 25 | Nylon | 0.45 | CS25N451 | 140.00 | CS25N452 | 635.00 | CS25N453 | 1,050.00 |
| 25 | PTFE | 0.20 | CS25T201 | 140.00 | CS25T202 | 635.00 | CS25T203 | 1,050.00 |
| 25 | PTFE | 0.45 | CS25T451 | 140.00 | CS25T452 | 635.00 | CS25T453 | 1,050.00 |
| 25 | PVDF | 0.45 | CS25P451 | 140.00 | CS25P452 | 635.00 | CS25P453 | 1,050.00 |
| 13 | Nylon | 0.20 | CS13N201 | \$120.00 | CS13N202 | \$535.00 | CS13N203 | \$920.00 |
| 13 | Nylon | 0.45 | CS13N451 | 120.00 | CS13N452 | 535.00 | CS13N453 | 920.00 |
| 13 | PTFE | 0.20 | CS13T201 | 120.00 | CS13T202 | 535.00 | CS13T203 | 920.00 |
| 13 | PTFE | 0.45 | CS13T451 | 120.00 | CS13T452 | 535.00 | CS13T453 | 920.00 |
| 13 | PVDF | 0.20 | CS13P201 | 120.00 | CS13P202 | 535.00 | CS13P203 | 920.00 |
| 13 | PVDF | 0.45 | CS13P451 | 120.00 | CS13P452 | 535.00 | CS13P453 | 920.00 |
| 4 | Nylon | 0.20 | CS04N201 | \$106.40 | CS04N202 | \$486.20 | CS04N203 | \$810.40 |
| 4 | Nylon | 0.45 | CS04N451 | 102.70 | CS04N452 | 459.30 | CS04N453 | 810.40 |
| 4 | PTFE | 0.20 | CS04T201 | 106.40 | CS04T202 | 486.20 | CS04T203 | 810.40 |
| 4 | PTFE | 0.45 | CS04T451 | 106.40 | CS04T452 | 486.20 | CS04T453 | 810.40 |

Pricing is subject to change without prior notice. For the most current pricing, contact the CSI Customer Service team at 1-800-267-8103 or email sales@chromspec.com