


Frustrated with leaky, hard to use press-fit connectors?



Now anyone in your lab can easily connect guard columns or retention gaps to GC columns by using Agilent Technologies' **Ultimate Union kits.**

- Ensure a reliable and high-quality connection—for leak-free results, time after time.
- Extend the life of your GC column—quickly and easily connect a guard column to your GC column.
- Reuse your column union for a dependable connection with the lowest dead volume union available, giving you confidence in your results.



Ultimate Union kits eliminate the “arts and crafts” involved in connecting columns with a simple-to-use, fast and reliable solution, helping your lab work smarter.



Agilent Technologies

Agilent Ultimate Union kits:

The smartest way to connect GC columns

- Easy to use
- Increased reliability
- Faster connections
- Longer column life

Ultimate Union Kit (G3182-61580) – Deactivated
Ultimate Union Kit (G3182-61581) – Non-deactivated

Each Agilent Ultimate Union kit contains:

- 1 union (your choice of deactivated or non-deactivated)
- 2 internal nuts
- 1 swaging nut
- 1 oven wall clip

Please remember to order ferrules suitable for your column diameter:

Silitite Ferrule Description	Quantity	Part Number
For 100-250 μm ID columns	10/pack	5188-5361
For 320 μm ID columns	10/pack	5188-5362
For 530 μm ID columns	10/pack	5188-5363



To order or learn more about Ultimate Union column connectors:
Contact your local Agilent representative or Agilent Authorized Distributor or visit www.agilent.com/chem/ultimateunion

This information is subject to change without notice.
© Agilent Technologies, Inc. 2005
Printed in U.S.A. October 31, 2005



Agilent Technologies