

Save up to  
**30%**  
 and get a FREE  
 inlet liner poster

## Inlet Liner Special Offer

Save 10%, 20%, or 30% off select inlet liners and receive a free technical liner poster.

Buy one 5-pack and receive a **10% discount** and a free technical liner poster.

Buy two\* 5-packs and receive a **20% discount** and a free technical liner poster.

Buy three\* 5-packs and receive a **30% discount** and a free technical liner poster.

Simply place your order referencing the catalog number from the following table, but hurry, offer expires September 30, 2008.

\*Discount applies to multiple packs of the same liner catalog number from the table below. Cannot be combined with other discounts or offers.

### For Agilent GCs

Description	ID* x OD x Length	qty.	cat.#
<b>Cyclo splitter® Liner for Agilent GCs</b>			
Cyclo splitter, 10% off	4.0mm x 6.3mm x 78.5mm	5-pk.	20707-409.1
Cyclo splitter, 20% off	4.0mm x 6.3mm x 78.5mm	2 x 5-pk.	20707-409.2
Cyclo splitter, 30% off	4.0mm x 6.3mm x 78.5mm	3 x 5-pk.	20707-409.3
<b>Cup Splitter Liner for Agilent GCs</b>			
Cup Splitter, 10% off	4.0mm x 6.3mm x 78.5mm	5-pk.	20710-409.1
Cup Splitter, 20% off	4.0mm x 6.3mm x 78.5mm	2 x 5-pk.	20710-409.2
Cup Splitter, 30% off	4.0mm x 6.3mm x 78.5mm	3 x 5-pk.	20710-409.3
<b>Drilled Uniliner® Liner for Agilent GCs (hole near bottom)</b>			
Drilled Uniliner (hole near bottom), 10% off	4.0mm x 6.3mm x 78.5mm	5-pk.	20771-409.1
Drilled Uniliner (hole near bottom), 20% off	4.0mm x 6.3mm x 78.5mm	2 x 5-pk.	20771-409.2
Drilled Uniliner (hole near bottom), 30% off	4.0mm x 6.3mm x 78.5mm	3 x 5-pk.	20771-409.3
<b>Siltek Drilled Uniliner (hole near bottom), 10% off</b>			
Siltek Drilled Uniliner (hole near bottom), 20% off	4.0mm x 6.3mm x 78.5mm	2 x 5-pk.	20771-214409.2
Siltek Drilled Uniliner (hole near bottom), 30% off	4.0mm x 6.3mm x 78.5mm	3 x 5-pk.	20771-214409.3
<b>Drilled Uniliner® Liner for Agilent GCs (hole near top)</b>			
Drilled Uniliner (hole near top), 10% off	4.0mm x 6.3mm x 78.5mm	5-pk.	21055-409.1
Drilled Uniliner (hole near top), 20% off	4.0mm x 6.3mm x 78.5mm	2 x 5-pk.	21055-409.2
Drilled Uniliner (hole near top), 30% off	4.0mm x 6.3mm x 78.5mm	3 x 5-pk.	21055-409.3
<b>Siltek Drilled Uniliner (hole near top), 10% off</b>			
Siltek Drilled Uniliner (hole near top), 20% off	4.0mm x 6.3mm x 78.5mm	2 x 5-pk.	21055-214409.2
Siltek Drilled Uniliner (hole near top), 30% off	4.0mm x 6.3mm x 78.5mm	3 x 5-pk.	21055-214409.3
<b>Split Liner w/ Wool for Agilent GCs</b>			
4mm Split w/ Wool, 10% off	4.0mm x 6.3mm x 78.5mm	5-pk.	20782-409.1
4mm Split w/ Wool, 20% off	4.0mm x 6.3mm x 78.5mm	2 x 5-pk.	20782-409.2
4mm Split w/ Wool, 30% off	4.0mm x 6.3mm x 78.5mm	3 x 5-pk.	20782-409.3
<b>Splitless Liner for Agilent GCs</b>			
4mm Splitless, 10% off	4.0mm x 6.5mm x 78.5mm	5-pk.	20773-409.1
4mm Splitless, 20% off	4.0mm x 6.5mm x 78.5mm	2 x 5-pk.	20773-409.2
4mm Splitless, 30% off	4.0mm x 6.5mm x 78.5mm	3 x 5-pk.	20773-409.3
<b>Gooseneck Splitless Liner for Agilent GCs</b>			
Gooseneck Splitless (4mm), 10% off	4.0mm x 6.5mm x 78.5mm	5-pk.	20799-409.1
Gooseneck Splitless (4mm), 20% off	4.0mm x 6.5mm x 78.5mm	2 x 5-pk.	20799-409.2
Gooseneck Splitless (4mm), 30% off	4.0mm x 6.5mm x 78.5mm	3 x 5-pk.	20799-409.3
<b>Siltek Gooseneck Splitless (4mm), 10% off</b>			
Siltek Gooseneck Splitless (4mm), 20% off	4.0mm x 6.5mm x 78.5mm	2 x 5-pk.	20799-214409.2
Siltek Gooseneck Splitless (4mm), 30% off	4.0mm x 6.5mm x 78.5mm	3 x 5-pk.	20799-214409.3
<b>Gooseneck Splitless Liner w/ Wool for Agilent GCs</b>			
Gooseneck Splitless (4mm) w/ Wool, 10% off	4.0mm x 6.5mm x 78.5mm	5-pk.	22406-409.1
Gooseneck Splitless (4mm) w/ Wool, 20% off	4.0mm x 6.5mm x 78.5mm	2 x 5-pk.	22406-409.2
Gooseneck Splitless (4mm) w/ Wool, 30% off	4.0mm x 6.5mm x 78.5mm	3 x 5-pk.	22406-409.3

\*Nominal ID at needle syringe expulsion point.



Use our "Liner Anatomy" poster to understand liner features and quickly determine which inlet liner is right for your application. Easy-to-read, informative comparisons of packing options, geometries, deactivations, and dimensions.

**FREE** with your inlet liner purchase from this flyer!

Save up to  
**30%**  
and get a **FREE**  
inlet liner poster

### For PerkinElmer GCs

Description	ID* x OD x Length	qty.	cat.#
<b>Cyclo-Uniliner® Liner for PerkinElmer GCs</b>			
Auto SYS Cyclo-Uniliner, 10% off	4.0mm x 6.2mm x 92.1mm	5-pk.	20840-409.1
Auto SYS Cyclo-Uniliner, 20% off	4.0mm x 6.2mm x 92.1mm	2 x 5-pk.	20840-409.2
Auto SYS Cyclo-Uniliner, 30% off	4.0mm x 6.2mm x 92.1mm	3 x 5-pk.	20840-409.3



<b>Split Precision™ Liner w/ Wool for PerkinElmer GCs</b>			
Auto SYS Split Precision Liner w/ Wool, 10% off	4.0mm x 6.2mm x 92.1mm	5-pk.	21027-409.1
Auto SYS Split Precision Liner w/ Wool, 20% off	4.0mm x 6.2mm x 92.1mm	2 x 5-pk.	21027-409.2
Auto SYS Split Precision Liner w/ Wool, 30% off	4.0mm x 6.2mm x 92.1mm	3 x 5-pk.	21027-409.3



<b>Split/Splitless Liner w/ Wool for PerkinElmer GCs</b>			
Auto SYS XL PSS Split/Splitless w/ Wool, 10% off	2.0mm x 4.0mm x 86.2mm	5-pk.	21718-409.1
Auto SYS XL PSS Split/Splitless w/ Wool, 20% off	2.0mm x 4.0mm x 86.2mm	2 x 5-pk.	21718-409.2
Auto SYS XL PSS Split/Splitless w/ Wool, 30% off	2.0mm x 4.0mm x 86.2mm	3 x 5-pk.	21718-409.3



### For Shimadzu GCs

Description	ID* x OD x Length	qty.	cat.#
<b>Split/Splitless Liner w/ Wool for Shimadzu GCs</b>			
Split/Splitless w/ Wool, 10% off	3.5mm x 5.0mm x 95mm	5-pk.	20956-409.1
Split/Splitless w/ Wool, 20% off	3.5mm x 5.0mm x 95mm	2 x 5-pk.	20956-409.2
Split/Splitless w/ Wool, 30% off	3.5mm x 5.0mm x 95mm	3 x 5-pk.	20956-409.3
<b>Siltek Split/Splitless w/ Wool</b>			
Siltek Split/Splitless w/ Wool, 10% off	3.5mm x 5.0mm x 95mm	5-pk.	20956-213409.1
Siltek Split/Splitless w/ Wool, 20% off	3.5mm x 5.0mm x 95mm	2 x 5-pk.	20956-213409.2
Siltek Split/Splitless w/ Wool, 30% off	3.5mm x 5.0mm x 95mm	3 x 5-pk.	20956-213409.3



<b>Uniliner® Liner w/ Wool for Shimadzu GCs</b>			
Uniliner w/ Wool, 10% off	3.5mm x 5.0mm x 95mm	5-pk.	21719-409.1
Uniliner w/ Wool, 20% off	3.5mm x 5.0mm x 95mm	2 x 5-pk.	21719-409.2
Uniliner w/ Wool, 30% off	3.5mm x 5.0mm x 95mm	3 x 5-pk.	21719-409.3



### For Thermo Scientific GCs

Description	ID* x OD x Length	qty.	cat.#
<b>Split/Splitless Siltek® Treated Metal Inlet Liner for Thermo Scientific GCs</b>			
Split/Splitless Siltek Metal Inlet Liners, 10% off	2.0mm x 2.75mm x 120mm	5-pk.	21083-409.1
Split/Splitless Siltek Metal Inlet Liners, 20% off	2.0mm x 2.75mm x 120mm	2 x 5-pk.	21083-409.2
Split/Splitless Siltek Metal Inlet Liners, 30% off	2.0mm x 2.75mm x 120mm	3 x 5-pk.	21083-409.3



<b>Split Liner for Thermo Scientific GCs</b>			
5mm Split, 10% off	5.0mm x 8.0mm x 105mm	5-pk.	20940-409.1
5mm Split, 20% off	5.0mm x 8.0mm x 105mm	2 x 5-pk.	20940-409.2
5mm Split, 30% off	5.0mm x 8.0mm x 105mm	3 x 5-pk.	20940-409.3



<b>Laminar Cup Splitter Liner for Thermo Scientific GCs</b>			
Laminar Cup Splitter, 10% off	4.0mm x 8.0mm x 105mm	5-pk.	20949-409.1
Laminar Cup Splitter, 20% off	4.0mm x 8.0mm x 105mm	2 x 5-pk.	20949-409.2
Laminar Cup Splitter, 30% off	4.0mm x 8.0mm x 105mm	3 x 5-pk.	20949-409.3

