

## CHROMSPEC™ CHROMATOGRAPHIC GRADE TUBING

The material is furnished with the ID free from oil, grease and foreign material. Careful quality control insures accuracy of size, soundness of material and uniform surface.

Product Number	Size	Description	Price
<b>Stainless Steel (Premium Quality, Chromatographic Grade)</b>			
C49300	1/16"	OD x 0.010" ID 316SS (Any length up to 300 ft.)	\$4.90 per ft.
C49305	1/16"	OD x 0.030" ID 316SS (Any length up to 300 ft.)	5.70 per ft.
C49306	1/16"	OD x 0.040" ID 316SS (Any length up to 300 ft.)	10.90 per ft.
C49310	1/8"	OD - 0.020 Wall 304SS (Up to 10 ft., straight lengths) (Min. 6 pcs.)	7.40 per ft.
C49315	1/8"	OD - 0.020 Wall 304SS 50 ft. continuous coils	263.30 per 50 ft.
C49320	3/16"	OD - 0.020 Wall 304SS (Up to 10 ft., straight lengths) (Min. 6 pcs.)	8.40 per ft.
C49325	1/4"	OD - 0.030 Wall 304SS (Up to 10 ft., straight lengths) (Min. 6 pcs.)	7.70 per ft.
<b>Nickel 200</b>			
C49330	1/16"	OD x 0.020" ID	8.10 per ft.
C49335	1/8"	OD Wall thickness 0.020"	9.70 per ft.
C49336	1/8"	OD Wall thickness 0.020"	350.50 per 50 ft.
C49340	1/4"	OD Wall thickness 0.020"	13.20 per ft.
<b>Copper (Refrigeration Grade Cleaned)</b>			
C49345	1/8"	OD Wall thickness 0.030"	68.50 per 50 ft.
C49350	3/16"	OD Wall thickness 0.030"	68.50 per 50 ft.
C49355	1/4"	OD Wall thickness 0.030"	68.50 per 50 ft.
<b>Aluminum</b>			
C49360	1/8"	OD Wall thickness 0.025"	90.20 per 50 ft.
C49370	1/4"	OD Wall thickness 0.032"	75.60 per 50 ft.
<b>PTFE Lined (Available only in 12 ft., 25 ft., and 50 ft. lengths)</b>			
C49395	1/8"	OD Stainless Steel Tubing PTFE Lined	17.00 per ft.
C49400	1/4"	OD Stainless Steel Tubing PTFE Lined	21.00 per ft.
<b>TFE/FEP</b>			
C49405	1/8"	OD x 0.084" ID Tubing FEP (Max. Temp. 204 °C)	210.30 per 100 ft.
C49410	1/8"	OD x 0.095" ID Tubing TFE (Max. Temp. 260 °C)	194.70 per 100 ft.
C49415	1/4"	OD x 3/16" ID Tubing FEP	126.90 per 25 ft.
C49420	1/4"	OD x 1/8" ID Tubing TFE	131.10 per 25 ft.

Other tubing materials and dimensions are available. Contact our helpful Technical Support Team for product numbers and availability.

Pricing is subject to change without prior notice. For the most current pricing, contact the CSI Customer Service team at 1-800-267-8103 or email [sales@chromspec.com](mailto:sales@chromspec.com)